

ZytoDot® 2C SPEC DDIT3 Break Apart Probe

Previously: ZytoDot 2C SPEC CHOP Break Apart Probe

Background


The ZytoDot® 2C SPEC DDIT3 Break Apart Probe is designed to detect translocations involving the chromosomal region 12q13.3 harboring the DDIT3 (C/EBP-homologous protein) gene (a.k.a. CHOP, GADD153) in formalin-fixed, paraffin-embedded tissue sections or cell samples. The DDIT3 gene encodes for a stress-induced dominant-negative inhibitor of the transcription factors C/EBP and LAP. DDIT3 is consistently rearranged in myxoid liposarcomas (MLS). The most frequent translocation involving the DDIT3 gene region is t(12;16)(q13.3;p11.2) and occurs in about 90% of patients with MLS. The rearrangement results in a fusion gene comprising the 5' part of the FUS (fused in sarcoma) gene, located in 16p11.2, and the complete coding region of the DDIT3 gene. The FUS-DDIT3 fusion protein acts as an abnormal transcription factor and development of myxoid liposarcomas is thus regarded as a consequence of deregulated FUS-DDIT3 target genes. Differential diagnosis of liposarcomas and accurate classification, the latter being especially important with regard to appropriate treatment and prognosis, are often problematic. Therefore, detection of DDIT3 rearrangements via ISH analysis is a valuable tool to confirm the histopathological diagnosis of myxoid liposarcoma.

References


- Aman P, et al. (1992) Genes Chromosomes Cancer 5: 278-85.
- Andersson M, et al. (2010) BMC Cancer 10: 249-58.
- Germano G, et al. (2010) Cancer Res 70: 2235-44.
- Meis-Kindblom JM, et al. (2001) Virchows Arch 439: 141-51.
- Panagopoulos I, et al. (1994) Cancer Res 54: 6500-3.
- Ron D & Habener JF (1992) Genes Dev 6: 439-53.

Probe Description

The ZytoDot® 2C DDIT3 Break Apart Probe is a mixture of a Digoxigenin-labeled and a Dinitrophenyl-labeled probe hybridizing to the 12q13.3-q14.1 band. The DNP-labeled probe hybridizes proximal to the DDIT3 gene and the DIG-labeled probe hybridizes distal to that gene.




Ideogram of chromosome 12 indicating the hybridization locations.




SPEC DDIT3 Probe map (not to scale).

Results

In an interphase nucleus of a normal cell lacking a translocation involving the 12q13.3-q14.1 band, using the ZytoDot® 2C CISH Implementation Kit, two red/green fusion signals are expected representing two normal (non-rearranged) 12q13.3-q14.1 loci. A signal pattern consisting of one red/green fusion signal, one red signal, and a separate green signal indicates one normal 12q13.3-q14.1 locus and one 12q13.3-q14.1 locus affected by a translocation or inversion.



SPEC DDIT3 Break Apart Probe hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.



Myxoid liposarcoma tissue section with translocation affecting the 12q13.3-q14.1 locus as indicated by one non-rearranged red/green fusion signal, one red signal, and one separate green signal indicating the translocation.

Prod. No.	Product	Label	Tests* (Volume)
C-3047-100	ZytoDot 2C SPEC DDIT3 Break Apart Probe CE IVD	Digoxigenin/DNP	10 (100 µl)
Related Products			
C-3044-10	ZytoDot 2C CISH Implementation Kit CE IVD Incl. Heat Pretreatment Solution EDTA, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 150 ml; 20x Wash Buffer TBS, 50 ml; Anti-DIG/DNP-Mix, 1 ml; HRP/AP-Polymer-Mix, 1 ml; AP-Red Solution A, 0.1 ml; AP-Red Solution B, 4 ml; HRP-Green Solution A, 0.2 ml; HRP-Green Solution B, 4 ml; Nuclear Blue Solution, 4 ml; Mounting Solution (alcoholic), 1 ml		10
WB-0009-500	Clear-it™ Stringency Buffer CE IVD		500 ml

* Using 10 µl probe solution per test. CE IVD only available in certain countries. All other countries research use only! Please contact your local dealer for more information.