

ZytoLight® SPEC SS18 Dual Color Break Apart Probe

Previously: ZytoLight SPEC SYT Dual Color Break Apart Probe

Background

The ZytoLight® SPEC SS18 Dual Color Break Apart Probe is designed to detect translocations involving the chromosomal region 18q11.2 harboring the SS18 (synovial sarcoma translocation, chromosome 18) gene (a.k.a. SYT).



Translocations involving the region 18q11.2 are found in over 90% of synovial sarcoma. Among soft tissue sarcomas, synovial sarcoma is one of the most common and classically occurs in the extremities of young adults with greater prevalence in males even though, the occurrence of synovial sarcoma has also been described in a wide variety of anatomical locations and in all ages. The most frequent translocation involving the SS18 gene region is t(X;18) (p11.23;q11.2) juxtaposing the SS18 gene in 18q11.2 either next to the SSX1 (synovial sarcoma, translocated to X chromosome) or the SSX2 gene, or very rarely to the SSX4 locus located in Xp11.23. Complex translocations involving other chromosomes are observed in less than 10% of synovial sarcomas. In combination with histopathological diagnosis, detection of SS18 rearrangements via FISH analysis is a valuable tool to confirm the diagnosis of synovial sarcoma.

References

- Amary MF, et al. (2007) *Mod Pathol* 20: 482-96.
 Clark J, et al. (1994) *Nat Genet* 7: 502-8.
 Ilmiawan MI, et al. (2012) *Med J Indones* 21: 196-202.
 Kawai A, et al. (1998) *N Engl J Med* 338: 153-60.
 Surace C, et al. (2004) *Lab Invest* 84: 1185-92.
 Torres L, et al. (2008) *Cancer Genet Cytogenet* 187: 45-9.


Probe Description

The SPEC SS18 Dual Color Break Apart Probe is a mixture of two direct labeled probes hybridizing to the 18q11.2 band. The orange fluorochrome direct labeled probe hybridizes distal to the SS18 gene, the green fluorochrome direct labeled probe hybridizes proximal to that gene.




Results

In an interphase nucleus lacking a translocation involving the 18q11.2 band two orange/green fusion signals are expected representing two normal (non-rearranged) 18q11.2 loci. A signal pattern consisting of one orange/green fusion signal, one orange signal, and a separate green signal indicates one normal 18q11.2 locus and one 18q11.2 locus affected by an 18q11.2 translocation.



SPEC SS18 Dual Color Break Apart Probe hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus.



Synovial sarcoma tissue section with translocation affecting the 18q11.2 locus as indicated by one non-rearranged orange/green fusion signal, one orange signal, and one separate green signal indicating the translocation.

Prod. No.	Product	Label	Tests* (Volume)
Z-2097-50	ZytoLight SPEC SS18 Dual Color Break Apart Probe CE IVD	●/●	5 (50 µl)
Related Products			
Z-2028-5	ZytoLight FISH-Tissue Implementation Kit CE IVD Incl. Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 150 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTest-Solution, 0.2 ml		5

* Using 10 µl probe solution per test. CE IVD only available in certain countries. All other countries research use only! Please contact your local dealer for more information.