



## INSTRUCTIONS:

### 1. Column preparation

Purge the pump with Binding Buffer removing all the air.

Connect the NeoTrap column to the pump by removing the end of the column and the top stop plug (save it for storage). Avoid introducing air in the column.

Wash the column with 10 column volumes of distilled water to eliminate the preservative.

### 2. Column equilibration

Equilibrate the column with 5 - 10 column volumes of binding buffer.

### 3. Sample application

It is recommendable to dilute the sample containing the immunoglobulins 1:1 with binding buffer to avoid ionic and pH changes.

All samples should be filtered through a 0.22  $\mu$ m filter in order to remove particles before applying it into the column.

### 4. Column washing

Wash with the binding buffer until the O.D. 280 nm reaches the baseline level again, normally 5-10 column volumes.

### 5. Purified immunoglobulin elution

Elution of the sample is performed at a pH below 3.0. In order to neutralize the eluted fractions, it is recommended the addition of 0.1 ml of neutralization buffer per ml of purified antibody.

### 6. Storage of the column

Put the top and bottom stop plugs in the column and keep at 2-8°C in 20% ethanol. **Do not freeze.**

For reference only

For Research Use Only. Not for Diagnostic or Therapeutic Use.