

First-in-Class Animal Origin Free Chemically Defined Media for the hMSC

Distributed by:
CliniSciences Group



Scope	<ul style="list-style-type: none"> • Mass production, research and drug development of MSC
Superiority	<ul style="list-style-type: none"> • Convenience for culture (ready to use) • No Coating material or supplement • Superior Proliferation
Application	<ul style="list-style-type: none"> • Adipose derived MSC • Bone marrow derived MSC • Umbilical cord derived MSC • Various origin MSC
competitiveness	<ul style="list-style-type: none"> • Stability of cell • Safety, homogeneity

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group


Proliferation

Test

- Test conducted to validate the media's proliferation rate

Methods

- Measure population doubling time (PDT)



Experimental Method: 1x10⁵ cells /25 T flask seeding, Sub-culture at 85-90% confluency conditions
Measuring and analyzing instrument: NC-250

The proliferation rate of MSCs cultured with CellCor CD AOF was faster and more stable across the field compared to fetal bovine serum, serum-free media, and our CellCor CD. In addition, the cumulative number of cells obtained was significantly higher for MSCs cultured with CellCor CD AOF (Donor N=3 or more).


Cell Viability

Test

- AO*DAPI (Acridine Orange*4'-6-Diamidino-2-phenylindole)

Methods

- Viability is measured by the ratio of total cells stained with AO and non-viable cells stained with DAPI.



Experimental Method: 1x10⁵ cells /25 T flask seeding, Sub-culture at 85-90% confluency conditions
Measuring and analyzing instrument: NC-250

All mesenchymal stem cells cultured with fetal bovine serum, serum-free medium, CellCor CD, and CellCor CD AOF were confirmed to have stable survival rates across all passages (Donor N=3 or more).

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Specific Marker

Test

- Purity (CD14, CD34, CD45)
- Identity (CD73, CD90, CD105)
- Immunogenicity (HLA-DR)

Methods

- Purity < 2%
- Identity > 95%
- Immunogenicity < 2%

Cell surface marker

	Surface marker	CellCor CD AOF	CellCor CD	FBS media	Commercial A
Purity (%)	CD14	0.5	0.2	0.7	10.3
	CD34	0.7	1.7	12.9	2.8
	CD45	0.6	0	0.7	0.9
Identity (%)	CD73	100	100	100	99.9
	CD90	100	100	99.9	99.9
	CD105	100	99.7	99.1	99.2
Immunogenicity (%)	HLA-DR	1.0	1.7	0.9	3.1

Experimental Method:

- 1x10⁵ cells/25 T flask seeding
- Analysis after 1 hour of staining with CD marker antibody
- Measurement and analysis instrument: Flow Cytometry


The surface markers were compared through flow cytometry after culturing ADSC with CellCor CD AOF, CellCor CD, fetal bovine serum, and serum-free medium. Specific surface markers that cultured with CellCor CD AOF were confirmed to be significantly maintained compared to other media (Donor N=3 or more, passage 10).

Differentiation

Test and Methods

- Adipogenesis : Oil red O Staining
- Chondrogenesis : Alcian blue Staining
- Osteogenesis : Alizarin red S Staining

Stain



Experimental Method:

- Adipogenesis: 2.1x10⁵ cells/12 well, Staining after 2 weeks of differentiation with StemPro™ Adipogenesis Differentiation Kit
- Chondrogenesis: 2.1x10⁴ cells/12 well, Staining after 2 weeks of differentiation with DMEM(High)+Ascorbic acid +Dexamethasone + TGF-β
- Osteogenesis: 2.1x10⁴ cells/12 well, Staining after 4 weeks of differentiation with StemPro™ Osteogenesis Differentiation Kit

Measuring and analysis instrument: NIKON TS2 microscope

When differentiated into Adipogenesis (adipocyte differentiation), Chondrogenesis (chondroblast differentiation), and Osteogenesis (osteoblast differentiation), respectively, after culturing ADSC with CellCor CD AOF, CellCor CD, fetal bovine serum, and serum-free medium, it was confirmed that ADSC differentiation ability was stably maintained compared to fetal bovine and serum-free medium (Donor N=3, passage 7).

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Senescence


Test

- Test that validates changes in cell senescence post sub-culture

Methods


- Identify dyed senescent cells to measure the amount of senescent cells

Cell size




- Experimental Method: 1×10^5 cells/25 T flask seeding, Sub-culture at 85-90% confluency conditions
- Measuring and analyzing instrument: NC-250

Stain



- Experimental Method: 3.8×10^4 cells/6-well seeding, After 3 days, staining with X-gal reagent, absorbance measured at 405nm (**p < 0.01)
- Measuring and analyzing instrument: Microscope, Microplate reader

Senescence gene expression



- Experimental Method: 1×10^5 cells/25 T flask seeding, After 3 days, harvest ADSC and qRT-PCR (***)
- Measuring and analyzing instrument: QuantStudio Real-Time PCR

ADSC cultured with CellCor CD AOF showed the lowest senescence (Donor N=3, passage 7)

Distributed by:

CliniSciences Group

Comparison of culture test with FBS and Commercial Medium

Animal/Serum Derivative Identification Test

Test



- Animal-derived serum (Neu5Gc) Identification
- Human-derived serum (BDKRB1) Identification

Methods

- Molecular Analysis of Neu5Gc Found in Mammals
- Identify BDKRB1 gene expression, a human-derived serum marker

Animal-derived serum (Neu5Gc) Identification

Human-derived serum (BDKRB1) Identification



Experimental Method: Measured after 1 hour of staining with Neu5Gc antibody on mesenchymal stem cells cultured in their respective culture media.

Measurement and analysis instrument: Flow cytometry

Experimental Method: 1×10^5 cells/25 T flask seeding/3 days later mesenchymal stem cell harvest, qRT-PCR (***) $p < 0.001$

Measurement and analyzer: QuantStudio Real-Time PCR

We checked animal-derived serum components and human-derived serum components in CellCor CD AOF and CellCor CD, fetal bovine serum, and serum-free media. We confirmed that CellCor CD AOF and CellCor CD media do not contain animal- or human-derived components (Donor N=3, passage 5).

Genetic Stability Comparison Test



Test

- Identify genetic mutations and abnormal DNA regulation during cell division

Methods

- Staining of mesenchymal cells to determine number and morphology

Karyotype



Experimental method: 1×10^5 cells/25 T flask seeding/staining of chromosomes of metaphase cells to check for abnormal mutations such as checking the shape, number, and arrangement of chromosomes.

Measuring and analyzing equipment: Request for Gendix analysis

We confirmed less difference in karyotype between lots in CellCor CD AOF and CellCor CD, Commercial A compared to MSC cultured in FBS media (Donor N=2, passage 7).

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Cell Self-Renewal


Test

- Identify Colony Forming assay

Methods

- Cationic dye, crystal violet, stains negatively charged DNA, proteins, etc. in the cell purple

Colony Forming assay



Experimental method: 3×10^2 cells/6 wells seeding/2 day medium change/2 weeks later staining with crystal violet

Measurement and analysis instrument: camera shooting

A colony forming assay was performed to determine self-renewal capacity after culturing ADSCs with CellCor CD AOF and CellCor CD, fetal bovine serum, and serum-free media. The highest self-renewal capacity was observed in ADSCs cultured with CellCor CD AOF.

Cell Tumorigenicity


Test

- Telomerase Reverse Transcriptase (TERT) Quantitative Test

Methods

- Sandwich ELISA assay for the quantitative measurement of Telomerase Reverse Transcriptase (450 nm)

TERT assay



Experimental Methods: Quantitative determination of Telomerase Reverse Transcriptase (450 nm) using culture medium after incubation of cells with each medium. (*p < 0.05, ***p < 0.001)

Measurement and analysis equipment: Microplate reader

A quantitative TERT assay was performed to determine the tumorigenicity of ADSCs cultured with CellCor CD AOF and CellCor CD, fetal bovine serum, and serum-free media. TERT expression was highest in ADSCs cultured with fetal bovine serum and lowest when cultured with CellCor CD AOF (N=3, passage 4).

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Genetic Stability

Test

- Identify STR profile

Methods

- Analyze whether STR differences occur due to medium.
- BIONICS outsourcing service

STR profile

Sample	Sample	% Match	Interpretation
1. A90 FBS P2	2. A90 FBS P7	100%	Related (Same donor)
3. A90 Commercial A P2	4. A90 StemPro P7	100%	Related (Same donor)
5. A90 CD P2	6. A90 CD P7	100%	Related (Same donor)
7. A90 AOF P2	8. A90 AOF P7	100%	Related (Same donor)

Experimental method: After genomic DNA extraction, STR locus amplification by PCR and STR pattern analysis

Measurement and analysis instrument: Applied biosystems simpliAmp PCR cycler, Applied biosystems 3730XL DNA analyzer

According to STR analysis, it was confirmed that the STR pattern was not altered even after culturing for 5 passages in all 4 types of media (Donor N=2, passages 2, 7).

	Sample AOF	Sample CD	Sample FBS	Sample Commercial A	Result
Marker	Allele	Allele	Allele	Allele	
D8S1179	11 12	11 12	11 12	11 12	Comparison
D21S11	30	30	30	30	
D7S820	10 12	10 12	10 12	10 12	
CSF1PO	9 10	9 10	9 10	9 10	
D3S1358	15 17	15 17	15 17	15 17	
TH01	6.3 7 8.3 9	6.3 7 8.3 9	6.3 7 8.3 9	6.3 7 8.3 9	
D13S317	11 12	11 12	11 12	11 12	
D16S539	9 10	9 10	9 10	9 10	
D2S1338	17 23	17 23	17 23	17 23	
D19S433	14	14	14	14	
Vwa	14 17	14 17	14 17	14 17	
TPOX	8 11	8 11	8 11	8 11	
D18S51	14	14	14	14	
AMEL	X Y	X Y	X Y	X Y	
D5S818	8 9 12 13	8 9 12 13	8 9 12 13	8 9 12 13	
FGA	22 23 24.2	22 23 24.2	22 23 24.2	22 23 24.2	

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Cell Culture Coating Verification Test


Test

- Verify Cell Proliferation.


Methods

- Verify the number of cells acquired

Cell Morphology



Acquired Cell Number



Experimental Method : 1x10⁵ cells/25 T flasks seeded and sub-cultured after 72 hours to determine cell number (***) p<0.001)

Measurement and analyzer : Nikon TS2 Microscope, Chemometec NC-250

We tested for cell culture with CellCor MSC CD AOF and Commercial A media. The Commercial A showed a difference in cell culture with and without coating, but CellCor MSC CD AOF worked well for cell culture regardless of coating (N=2, passage 4).

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Gene expression analysis and functional validation


Test

- Analysis of senescence related gene expression

Methods

- Identify of expression of senescence-related genes IGFBP5 (Insulin-like growth factor binding protein 5) and Endoglin

Senescence-related gene expression



Experimental Method: 1×10^5 cells/25 T flasks seeding/3 days later mesenchymal stem cell harvest, qRT-PCR confirmation, (*p < 0.01, ***p < 0.001)

Measurement and analyzer: QuantStudio Real-Time PCR

Senescence-related gene expression was highest in culture with fetal bovine serum medium and it showed lower expression in mesenchymal stem cells cultured with CellCor CD AOF (Donor N=3, passage 7).


Test

- Analysis of apoptosis related gene expression.

Methods

- Identify of expression of apoptosis-related gene G0G1

Apoptosis-related gene expression



Experimental Method: 1×10^5 cells/25 T flasks seeding/3 days later mesenchymal stem cell harvest, qRT-PCR confirmation, (***p < 0.001)

Measurement and analyzer: QuantStudio Real-Time PCR

Apoptosis-related gene expression was highest in culture with fetal bovine serum medium and it showed lower expression in mesenchymal stem cells cultured with CellCor CD AOF (Donor N=3, passage 7).

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Gene expression analysis and functional validation


Test

- Analyze immune response-related genes

Methods

- Identify expression of immune response-related genes CXCL12 (C-X-C motif chemokine 12), PTX3 (Pentraxin 3)

Immune response related gene expression



Experimental Method: 1×10^5 cells/25 T flasks seeding/3 days later mesenchymal stem cell harvest, qRT-PCR confirmation, (*p < 0.05, ***p < 0.001)
Measurement and analyzer: QuantStudio Real-Time PCR

Immune response-related gene expression was highest in culture with fetal bovine serum medium and it showed lower expression in mesenchymal stem cells cultured with CellCor CD AOF (Donor N=3, passage 7).


Test

- Analyze inflammatory response-related genes

Methods

- Identify expression of inflammatory response-related genes TNFRSF11B (tumor necrosis factor receptor superfamily member 11b), CSF1 (colony-stimulating factor 1).

Inflammatory response related gene expression



Experimental Method: 1×10^5 cells/25 T flasks seeding/3 days later mesenchymal stem cell harvest, qRT-PCR confirmation, (*p < 0.05, ***p < 0.001)
Measurement and analyzer: QuantStudio Real-Time PCR

Inflammatory response-related gene expression was highest in culture with fetal bovine serum medium and it showed lower expression in mesenchymal stem cells cultured with CellCor CD AOF (Donor N=3, passage 7).

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Protein expression analysis and functional validation


Test

- Apoptosis-related protein analysis

Methods

- Identify expression of apoptosis-related proteins filamin A and STAT1 (signal transduction and activator of transcription 1)

Apoptosis-related protein analysis



Experimental Method: 1×10^5 cells/25 T flasks seeded/3 days later, mesenchymal stem cells harvested, protein expression confirmed by Western blot.

Measurement and analyzer : Image J, iBright CL1000

Apoptosis-related protein expression was higher in ADSCs cultured in fetal bovine serum medium and lower in ADSCs cultured in CellCor CD AOF. (Donor N=2, passage 7).


Test

- Inflammatory response-related protein analysis

Methods

- Identify the expression of PTX3 (Pentraxin 3), an inflammatory response-related protein, and Optineurin

Inflammatory response-related protein analysis



Experimental Method: 1×10^5 cells/25 T flasks seeded/3 days later, mesenchymal stem cells harvested, protein expression confirmed by Western blot.

Measurement and analyzer : Image J, iBright CL1000

Inflammatory response-related protein expression was higher in ADSCs cultured with fetal bovine serum medium and lower in ADSCs cultured with CellCor CD AOF. (Donor N=2, passage 7)

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Protein expression analysis and functional validation


Test

- Oncogene candidate-related protein analysis

Methods

- Identify expression of oncogene candidate-related protein PDLIM5 (PDZ and LIM domain protein 5)

Oncogene candidate-related protein analysis



Experimental Method: 1×10^5 cells/25 T flasks seeded/3 days later, mesenchymal stem cells harvested, protein expression confirmed by Western blot.
Measurement and analyzer : Image J, iBright CL1000

Oncogene candidate-related protein expression was higher in ADSCs cultured in fetal bovine serum medium and lower in ADSCs cultured in CellCor CD AOF (Donor N=2, passage 7).


Test

- Cancer overexpression-related protein analysis

Methods

- Identify expression of cancer overexpression-related protein MYOF (Myoferlin)

Cancer overexpression-related protein analysis



Experimental Method: 1×10^5 cells/25 T flasks seeded/3 days later, mesenchymal stem cells harvested, protein expression confirmed by Western blot.
Measurement and analyzer : Image J, iBright CL1000

Cancer overexpression-related protein expression was higher in ADSCs cultured in fetal bovine serum medium and lower in ADSCs cultured in CellCor CD AOF (Donor N=2, passage 7).

Comparison of culture test with FBS and Commercial Medium

Distributed by:
CliniSciences Group

Protein expression analysis and functional validation


Test

- DNA repair-related protein analysis

Methods

- Identify expression of DNA repair-related protein MCM (minichromosome maintenance) 2, 3, 7

DNA repair-related protein analysis



Experimental Method: 1×10^5 cells/25 T flasks seeded/3 days later, mesenchymal stem cells harvested, protein expression confirmed by Western blot.

Measurement and analyzer : Image J, iBright CL1000

DNA repair-related protein expression was higher in ADSCs cultured in CellCor CD AOF and lower in ADSCs cultured in fetal bovine serum medium (Donor N=2, passage 7).


Test

- Stabilization of mRNA-related protein analysis

Methods

- Identify expression of stabilization of mRNA-related protein SNRNP70 (Small nuclear ribonucleoprotein 70 kDa)

Stabilization of mRNA-related protein analysis



Experimental Method: 1×10^5 cells/25 T flasks seeded/3 days later, mesenchymal stem cells harvested, protein expression confirmed by Western blot.

Measurement and analyzer : Image J, iBright CL1000

Stabilization of mRNA-related protein expression was higher in ADSCs cultured in CellCor CD AOF and lower in ADSCs cultured in fetal bovine serum medium (Donor N=2, passage 7).

Comparison of culture test with FBS and Commercial Medium

Distributed by:

CliniSciences Group

Protein expression analysis and functional validation



Test

- Malignancy-related protein analysis

Methods

- Identify expression of malignancy-related protein IMP3 (Insulin-like growth factor II mRNA binding protein 3)

Malignancy-related protein analysis



Experimental Method: 1×10^5 cells/25 T flasks seeded/3 days later, mesenchymal stem cells harvested, protein expression confirmed by Western blot.

Measurement and analyzer : Image J, iBright CL1000

Malignancy-related protein expression was higher in ADSCs cultured in fetal bovine serum medium and lower in ADSCs cultured in CellCor CD AOF (Donor N=2, passage 7).

Comparison of culture test with Commercial CD medium

Distributed by:
CliniSciences Group

Proliferation


Test

- Test conducted to validate the media's proliferation rate


Methods

- Measure population doubling time (PDT)

PDT



ACN



Experimental Method: 1x10⁵ cells /25 T flask seeding, Sub-culture at 85-90% confluency conditions
Measuring and analyzing instrument: NC-250

The proliferation rate of MSCs cultured with CellCor CD AOF was faster and more stable across the field compared to Commercial B, and our CellCor CD. In addition, the cumulative number of cells obtained was significantly higher for MSCs cultured with CellCor CD AOF.

Viability


Test

- AO*DAPI (Acridine Orange*4'-6-Diamidino-2-phenylindole)

Methods

- Viability is measured by the ratio of total cells stained with AO and non-viable cells stained with DAPI.

Viability



Experimental Method: 1x10⁵ cells /25 T flask seeding, Sub-culture at 85-90% confluency conditions
Measuring and analyzing instrument: NC-250

All mesenchymal stem cells cultured with CellCor CD AOF, CellCor CD and Commercial B were confirmed to have stable survival rates across all passages.

Comparison of culture test with Commercial CD medium

Distributed by:
CliniSciences Group

Specific Marker

Test

- Purity (CD14, CD34, CD45)
- Identity (CD73, CD90, CD105)
- Immunogenicity (HLA-DR)

Methods

- Purity < 2%
- Identity > 95%
- Immunogenicity < 2%

Cell surface marker

	Surface marker	CellCor CD AOF	CellCor CD	Commercial B
Purity (%)	CD14	0.5	0.01	0.05
	CD34	0.7	3.83	15.77
	CD45	0.6	0.01	0.08
Identity (%)	CD73	100	99.88	99.55
	CD90	100	99.69	98.30
	CD105	100	98.75	99.33
Immunogenicity (%)	HLA-DR	1.0	0.01	0.02

Experimental Method: 1x10⁵ cells/25 T flask seeding, Analysis after 1 hour staining with CD marker antibody
Measuring and analyzing instrument: Flow Cytometry


After culturing ADSC with CellCor CD AOF, CellCor CD and Commercial B, the surface markers were compared through flow cytometry. Specific surface markers that cultured with CellCor CD AOF were confirmed to be significantly maintained compared to other media (Passage 10).

Differentiation

Test and Methods

- Adipogenesis : Oil red O Staining
- Chondrogenesis : Alcian blue Staining
- Osteogenesis : Alizarin red S Staining

Stain



Experimental Method:

- Adipogenesis: 2.1x10⁵ cells/12 well, Staining after 2 weeks of differentiation with StemPro™ Adipogenesis Differentiation Kit
 - Chondrogenesis: 2.1x10⁴ cells/12 well, Staining after 2 weeks of differentiation with DMEM(High)+Ascorbic acid +Dexamethasone + TGF-β
 - Osteogenesis: 2.1x10⁴ cells/12 well, Staining after 4 weeks of differentiation with StemPro™ Osteogenesis Differentiation Kit
- Measuring and analyzing instrument: NIKON TS2 microscope

When differentiated into Adipogenesis (adipocyte differentiation), Chondrogenesis (chondroblast differentiation), and Osteogenesis (osteoblast differentiation), respectively, after culturing ADSC with CellCor CD AOF, CellCor CD and Commercial B. It was confirmed that ADSC differentiation ability was stably maintained compared to Commercial B (Passage 10).

Comparison of culture test with Commercial CD medium

Distributed by:
CliniSciences Group

Senescence



Test

- Test that validates changes in cell senescence post sub-culture

Methods

- Identify dyed senescent cells to measure the amount of senescent cells

Stain



Experimental Method: 3.8×10^4 cells/6-well seeding, After 3 days, staining with X-gal reagent
Measuring and analyzing instrument: Microscope

ADSC cultured with CellCor CD AOF showed the lowest senescence (Passage 6).

Human Serum Derivative Identification Test


Test

- Human-derived serum (BDKRB1) Identification

Methods

- Identify BDKRB1 gene expression, a human-derived serum marker

Human-derived serum (BDKRB1) Identification



Experimental Method: 1×10^5 cells/25 T flask seeding/3 days later mesenchymal stem cell harvest, qRT-PCR (***) $p < 0.001$

Measurement and analyzer: QuantStudio Real-Time PCR

We checked human-derived serum components in CellCor CD AOF, CellCor CD, Commercial A and Commercial B media. We confirmed that CellCor CD AOF, CellCor CD and Commercial B media do not contain animal-derived or human-derived components (Passage 6).

Comparison of culture test with Commercial CD medium

Distributed by:
CliniSciences Group

Cell Culture Coating Verification Test



Test

- Verify of cell attachment depending on coating material.

Methods

- After cell culturing 72 hours, analysis of the adhesion rate depending on coating materials.

Cell attachment



Experimental Method : 1x10⁵ cells/25 T flasks seeded and sub-cultured after 72 hours to determine cell number (***) p<0.001
Measurement and analyzer : NC-250

While there was no difference in cell culture between CellCor CD AOF and CellCor CD depending on the presence or absence of coating on the culture dish, we found a difference in cell culture of Commercial B on the presence or absence of coating (Passage 10).

Self-renewal


Test

- Identify Colony Forming assay

Methods

- Cationic dye, crystal violet, stains negatively charged DNA, proteins, etc. in the cell purple

Colony Forming assay



Experimental method: 3x10² cells/6 wells seeding/2 day medium change/2 weeks later staining with crystal violet
Measurement and analysis instrument: camera shooting

After culturing ADSCs with CellCor CD AOF, CellCor CD and Commercial B medium, a colony forming assay was performed to determine self-renewal capacity. The highest self-renewal capacity was observed in ADSCs cultured with CellCor CD AOF (Passage 2).

Comparison of culture test with various tissue origin

Distributed by:
CliniSciences Group

Cell morphology

Test


- Identify stem cell shapes by tissue origin

Methods

- Verify morphology after stem cell culture by tissue origin

Identify stem cell shapes by tissue origin

	Stem cell	Tissue
UCSC	Umbilical Cord-derived Stem Cell (UCSC)	Umbilical Cord
BMSC	Bone Marrow-derived Stem Cell (BMSC)	Bone Marrow



Experimental Method: 1x10⁵ cells/25 T flask seeding, cell shape imaging (4X) under 85-90% confluency conditions
Measurement and analysis instrument: Microscope

We checked stem cell shape by tissue origin in CellCor CD AOF and CellCor CD MSC. We confirmed that stem cells cultured in CellCor CD AOF and CellCor CD MSC were smaller in size and maintained a more spindle-like shape compared to stem cells cultured in fetal bovine serum medium (passage 5).

Cell proliferation and viability assays



Test

- ACN (Accumulated Cell Number)
- AO*DAPI (Acridine Orange*4'-6-Diamidino-2-phenylindole)



Methods

- (cumulative cell count/cells seeded) x number of cells obtained
- Measures viability as the ratio of total cells stained with AO to non-viable cells stained with DAPI

ACN



Viability



Experimental Method: 1x10⁵ cells/25 T flask seeding, sub-culture at 85-90% confluency conditions
Measuring and analyzing instrument: NC250

CellCor CD AOF and CellCor CD MSC were observed to have a higher proliferation rate compared to fetal bovine serum medium. The cell viability was stable across the field (UCSC, BMSC Donor N=2).

Comparison of culture test with various tissue origin

Distributed by:
CliniSciences Group

Specific Marker Expression

Test

- Purity (CD14, CD34, CD45)
- Identity (CD73, CD90, CD105)
- Immunogenicity (HLA-DR)

Methods

- Purity < 2%
- Identity > 95%
- Immunogenicity < 2%

Cell surface marker

UCSC	Surface marker	CellCor CD AOF	CellCor CD	FBS media	BMSC	Surface marker	CellCor CD AOF	CellCor CD	FBS media
Purity (%)	CD14	0.47	0.06	0.1	Purity (%)	CD14	0.35	0.02	0.28
	CD34	0.86	0.23	0.27		CD34	1.25	0.19	0.2
	CD45	0.25	0.12	0.1		CD45	0.24	0.02	0.24
Identity (%)	CD73	99.99	99.82	99.71	Identity (%)	CD73	99.98	99.9	98.64
	CD90	99.99	99.79	99.67		CD90	99.99	99.98	97.89
	CD105	99.77	99.42	99.73		CD105	99.93	99.74	97.21
Immunogenicity (%)	HLA-DR	0.14	0.52	0.15	Immunogenicity (%)	HLA-DR	0.15	0	29

Experimental Method: 1x10⁵ cells/25 T flask seeding, Analysis after 1 hour staining with CD marker antibody
Measuring and analyzing instrument: Flow Cytometry

After culturing UCSC and BMSC with CellCor CD AOF, CellCor CD and fetal bovine serum, the surface markers were compared through flow cytometry. Specific surface markers that cultured with CellCor CD AOF were confirmed to be significantly maintained compared to other media (UCSC, BMSC Donor N=2, Passage 7)

CliniSciences Group

Austria

Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>



Belgium

Company: CliniSciences S.R.L.
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>



Denmark

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>



Finland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>



France

Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre - France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>



Germany

Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>



Iceland

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>



Ireland

Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>



Italy

Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>



Netherlands

Company: CliniSciences B.V.
Address: Gaetano Martinolaan 85, 6229
GS Maastricht - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>



Norway

Company: CliniSciences AS
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, Netherlands
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>



Poland

Company: CliniSciences sp.Z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>



Portugal

Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>



Spain

Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 91 269 40 65
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>



Sweden

Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: sverige@clinisciences.com
Web: <https://www.clinisciences.com>



Switzerland

Company: CliniSciences AG
Address: Fracht Ost Flughafen
Kloten CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>



UK

Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>



USA

Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>

