



**Distributed by:**

**CliniSciences Group**

# CRISPR-Cas9

TOOLS & SERVICES FOR GENE EDITING


Screening | Knock-out | Knock-in | Cell Line Development



# Simple and Effective Genome Editing

CRISPR-Cas9 is a technique used for genome editing that is adapted from bacterial antiviral immune mechanisms. Bacteria capture and store DNA fragments from invading viruses within a region of their genome, and these CRISPR (Clustered Regularly Interspaced Short Palindromic Repeats) guide sequences help detect and protect the bacteria from future infections. When the CRISPR guide sequences detect an invading virus or DNA whose sequence is complementary to the CRISPR guide, the Cas9 (CRISPR-associated protein 9) nuclease is recruited to specifically cleave the invading DNA, resulting in its degradation.

This CRISPR-Cas9 system has been modified for use in mammalian cells. By introducing a guide sequence (sgRNA) specific for our gene of interest, we can either knock-out specific genes through introducing frame shift mutations via Non-Homologous End Joining (NHEJ), or generate knock-in mutations by additionally providing a template for Homologous Recombination (HR).



## Our Advantages



### Conducted In-House

- All CRISPR products and services are produced and performed in the USA at our San Diego, CA laboratory
- Get customized, personal support directly from our CRISPR experts



### Committed to Excellence

- ISO 9001:2015-certified Quality Management System



### Customized For Your Research Needs

- Screening & Profiling: >200 optimized cell lines and cell-based assays
- Cell line development: choose from >30 cell types and >70 reporter systems
- Ready-to-use Lentiviruses: integrating and non-integrating options



### Multiple CRISPR Editing Tools & Applications

- CRISPR knock-out
- CRISPR knock-in
- CRISPR activation
- Stable cell lines, cell pools, lentiviruses, & plasmids
- CRISPR screens

Distributed by:


**1** **iniSciences Group**

# Knock-in Cell Lines

CRISPR-Cas9 can be used to introduce specific changes within the sequence of a gene, resulting in targeted protein alteration, a process known as gene knock-in. The knock-in can either be a single nucleotide substitution or an extended sequence encoding a full protein. This is an effective approach to study known protein mutations, to screen for mutations that affect protein function, or add tags to endogenous proteins in cell lines, in addition to other applications.

With BPS Bioscience's custom cell line development services, our team of highly experienced scientists can generate custom knock-in cell lines in more than 30 different cell types (or your preferred cell types) using CRISPR-Cas9 licensed technology, targeting your gene(s) of interest. The development process is comprised of distinct milestones where data is provided after each milestone completion. Each project can be fully customizable for your desired deliverables.

## A Milestone-Measured Process Ensures Success



BPS Bioscience will design and construct the sgRNA and HDR template according to your experimental needs.



The cells will be transduced with the Cas9, sgRNA, and HDR template, followed by genome editing evaluation. Single clones will be selected and expanded.




The knock-in mutations will be confirmed by genomic sequencing, and positive clones will be expanded for further confirmation.

## Consider the Advantages and Limitations of CRISPR-Cas9 Technology

### Advantages

- Precise gene editing
- Relatively high efficiency
- Permanent gene editing to generate stable knock-ins
- Relatively simple compared to other gene editing techniques
- Can target multiple genes simultaneously
- Effective across many cell types
- Potential for use in screening



### Limitations

- Potential for off-target effects
- Variable effectiveness of sgRNAs
- Target selection may be limited by PAM (protospacer adjacent motif) sequence
- Potential for low homology directed repair efficiency in your target cell
- Large proteins may be difficult to introduce

**Distributed by:** The milestone-measured process ensures that any potential limitations are overcome before advancing in the project.

**CliniSciences Group**

# Knock-out Cell Lines & Cell Pools

CRISPR-Cas9 is an ideal system for targeted gene knock-out. The achieved result can be a crude cell pool, in which the cell population is heterogenous with differing degrees of gene knock-out, or a stable cell line derived from a single clone. Stable cell lines are best for long-term or complex studies to provide experimental consistency over time. Cell pools are useful for lower cost, initial testing.


BPS Bioscience provides ready-to-use Cas9-expressing cell lines and cell pools that allow you to perform your own CRISPR knock-out assays. Alternatively, we can generate knock-outs for you with our milestone-measured process.

## Project Milestones for Knock-out Results



1  
Molecular  
Biology

We will synthesize up to 5 sgRNA sequences and clone into a CRISPR expression vector.




2  
CRISPR  
Transfection

Depending on the cell type, cells can be transduced via electroporation, liposome-based transfection, or viral infection.




3  
Limiting  
Dilution

Based on the results of the initial pool testing, the cell pool will be clonally diluted and single cell-derived clones will be expanded.



4  
Confirmation of  
Expression

The expression level of the target protein will be analyzed via Western blot or flow cytometry.




5  
Confirmation &  
Delivery

Gene knock-out will be confirmed through genomic sequencing. Confirmed clones will be expanded, frozen, and tested for mycoplasma contamination.

## Why knock-out genes in cells?

- Study transcriptional regulation by knocking out transcription factors, repressors, or epigenetic enzymes
- Generate disease models
- Generate knock-out libraries for screening
- Generate unique immune phenotypes, such as MHC-deficient cells



- Identify potential cancer targets
- Identify critical factors in signaling pathways
- Create systems for cell & gene therapy
- Optimize CAR-T cell function and limit toxicity
- Identify factors for viral entry
- And many more...

**Distributed by:**


3 **CliniSciences Group**

# CRISPR Activation


The CRISPR (SAM) system (CRISPR-based Synergistic Activation Mediator) is a combination of dCas9 and other molecular biology tools designed to activate the transcription of any endogenous gene of interest.

The system comprises 3 components that form a DNA-binding complex upon transfection into the cells. The first component is dCas9 (dead Cas9 with a disabled endonuclease activity) fused to transcriptional activator VP64, typically composed of four tandem copies of VP16 (Herpes Simplex Viral Protein 16, amino acids 437-447). The other two components exploit the unique MS2 bacteriophage protein/RNA interaction system in which the coat protein of the bacteriophage binds tightly and specifically to a distinct 19-nucleotide RNA aptamer. In the second component of SAM, MS2 aptamers forming a characteristic stem loop structure are added to the sgRNA. The sgRNA-MS2 component forms a complex with dCas9 and directs it to the target DNA sequence next to the promoter region of

the gene of interest. The sgRNA-MS2 aptamer recruits the third SAM component consisting of transcriptional activators P65 (Nuclear Factor NF- $\kappa$ B p65) and HSF1 (Heat Shock Factor 1) fused with an MS2-tag corresponding to the minimal aptamer-binding peptide of the MS2 coat protein. Once captured in the assembled complex at the gene promoter, P65 and HSF1 synergize with VP64 to robustly activate transcription of the downstream target gene, as much as a hundred-fold depending on the gene. Theoretically, the SAM system can be used to target one or several gene promoters in the same cell.




## Accelerate your discoveries with our CRISPRa solutions



### Cell Lines

- Stably expressing dCas9-VP64 and MS2-P65-HSF1
- Convenient and cost-effective platform for setting up your own activating experiments or screens



### Lentiviruses

- Integrating lentiviruses containing dCas9-VP64 and MS2-P65-HSF1
- Integrating sgRNA-MS2 lentiviruses, containing 4 validated sgRNA targeting your gene of interest
- Ready for transduction into almost all types of mammalian cells, including primary and non-dividing cells



### Plasmids

- sgRNA-MS2 plasmids can be transfected or electroporated into your cells for transient or stable (following drug selection) activation
- Virus-free option

**Distributed by:**

**CliniSciences Group**

AAVs	Catalog#	Cell Lines	Catalog#
AAV-DJ SaCas9	78478	CRISPRa (SAM) HeLa Cell Line	78193
AAV1 SaCas9	78479	CRISPRa (SAM) HepG2 Cell Line	78194
AAV2 SaCas9	78480	CRISPRa (SAM) Jurkat Cell Line	78080
AAV3 SaCas9	78481	CRISPRa (SAM) MCF7 Cell Line	78522
AAV5 SaCas9	78483	CRISPRa (SAM) MDA-MB-231 Cell Line	78521
AAV6 SaCas9	78484	FcGR1a (CD64) Knockout THP-1 Cell Line	82191
AAV8 SaCas9	78485	FCGR2A (CD32A) Knockout Jurkat Cell Line	78549
AAV9 SaCas9	78486	Firefly Luciferase CD19 Knockout NALM6 Cell Line	82168
		Firefly Luciferase CD19 Knockout Raji Cell Line	82167
		PARG Knockout HeLa Cell Line	82171
		PARP1 Knockout HeLa Cell Line	82169
		TCR Knockout Jurkat Cell Line	78539
		TCR Knockout NFAT-Luciferase Reporter Jurkat Cell Line	78556
		TCR/B2M Knockout NFAT Luciferase Reporter Jurkat Cell Line	78557
Cell Lines	Catalog#	Lentiviruses	Catalog#
B2M Knockout iPS Cell Line	82161	B2M (Human) CRISPR/Cas9 Lentivirus (Integrating)	78340
B2M Knockout Jurkat Cell Line	78342	B2M (Human) CRISPR/Cas9 Lentivirus (Non-Integrating)	78341
B2M Knockout NFAT Luciferase Reporter Jurkat Cell Line	78363	BCMA CRISPR/Cas9 Lentivirus (Integrating)	78893
B2M Knockout THP-1 Cell Line	78389	BCMA CRISPR/Cas9 Lentivirus (Non-Integrating)	78894
B2M/CIITA Double Knockout THP-1 Cell Line	78391	Cas9 Lentivirus (Hygromycin Selection)	78067
Cas9 Expressing A549 Cell Pool	78072	Cas9 Lentivirus (Neomycin Selection)	78432
Cas9 Expressing Daudi Cell Pool	78089	Cas9 Lentivirus (Puromycin Selection)	78066
Cas9 Expressing HCT116 Cell Pool	78073	CBL-B (Human) CRISPR/Cas9 Lentivirus (Integrating)	78343
Cas9 Expressing iPS Cell Pool	78578	CBL-B (Human) CRISPR/Cas9 Lentivirus (Non-Integrating)	78344
Cas9 Expressing Jurkat Cell Pool	78070	CD47 CRISPR/Cas9 Lentivirus (Integrating)	78056
Cas9 Expressing MDA-MB-231 Cell Pool	78069	CD47 CRISPR/Cas9 Lentivirus (Non-Integrating)	78063
Cas9 Expressing Raji Cell Pool	78071	CD5 (Human) CRISPR/Cas9 Lentivirus (Integrating)	78119
Cas9 Inducible (Tet-On) iPS Cell Pool	78845	CD5 (Human) CRISPR/Cas9 Lentivirus (Non-Integrating)	78198
Cas9-Expressing A549 Cell Line (High Expression or Low Expression)	78134	CIITA (Human) CRISPR/Cas9 Lentivirus (Integrating)	78435
Cas9-Expressing Daudi Cell Line	78157	CIITA (Human) CRISPR/Cas9 Lentivirus (Non-integrating)	78434
Cas9-Expressing HCT116 Cell Line (High or Low Expression)	78135	CRBN CRISPR/Cas9 Lentivirus (Integrating)	78517
Cas9-Expressing HEK293 Cell Line	78166	CRBN CRISPR/Cas9 Lentivirus (Non-Integrating)	78518
Cas9-Expressing HeLa Cell Pool	78161	CRISPR/Cas9 Kinase Knockout Lentivirus Library (Array Format)	78487
Cas9-Expressing Jurkat Cell Line (High or Low Expression)	78136	CTLA4 CRISPR/Cas9 Lentivirus (Integrating)	78054
Cas9-Expressing MCF7 Cell Pool	78179	CTLA4 CRISPR/Cas9 Lentivirus (Non-Integrating)	78061
Cas9-Expressing MDA-MB-231 Cell Line (High or Low Expression)	78150	FCGR2A CRISPR/Cas9 Lentivirus (Integrating)	78537
Cas9-Expressing Neuro2A Cell Line (High or Low Expression)	78137	FCGR2A CRISPR/Cas9 Lentivirus (Non-Integrating)	78538
Cas9-Expressing Neuro2A Cell Pool	78087	Kinase (Human) CRISPR/Cas9 Lentivirus (Integrating)	78488
Cas9-Expressing Raji Cell Line	78156	LAG3 CRISPR/Cas9 Lentivirus (Integrating)	78053
Cas9/GFP Safe-Harbor HEK293 Cell Line	78582		
CD19 Knockout Raji Cell Line	82166		
CD36 Knockout THP-1 Cell Line	82190		
CD8+ TCR Knockout NFAT-Luciferase Reporter Jurkat Cell Line	78757		
CIITA Knockout THP-1 Cell Line	78390		
CRISPRa (SAM) HEK293 Cell Line	78192		

**Distributed by:**

**CliniSciences Group**

Lentiviruses	Catalog#
LAG3 CRISPR/Cas9 Lentivirus (Non-Integrating)	78060
PD-1 (Human) sgRNA-MS2 Lentivirus (Integrating)	78190
PD-1 CRISPR/Cas9 Lentivirus (Integrating)	78052
PD-1 CRISPR/Cas9 Lentivirus (Non-Integrating)	78059
PD-L1 CRISPR/Cas9 Lentivirus (Integrating)	78057
PD-L1 CRISPR/Cas9 Lentivirus (Non-Integrating)	78064
TCR CRISPR/Cas9 Lentivirus (Integrating)	78055
TCR CRISPR/Cas9 Lentivirus (Non-Integrating)	78062
TGFBR2 CRISPR/Cas9 Lentivirus (Integrating)	78535
TGFBR2 CRISPR/Cas9 Lentivirus (Non-Integrating)	78536
TIGIT CRISPR/Cas9 Lentivirus (Integrating)	78058
TIGIT CRISPR/Cas9 Lentivirus (Non-Integrating)	78065

Proteins	Catalog#
Cas12, GST-tag (Lachnospiraceae)	100740
Cas12, His-tag (Lachnospiraceae)	100741
Cas12a, His-Tag (Acidaminococcus sp.)	101627
Cas12b (A. acidiphilus)	101626
Cas13a (L. buccalis)	101629
Cas13a, His-Tag (L. wadei)	101630
Cas13b, His-Tag (Prevotella sp)	101631
Cas9 (D10A), NLS, His-Tag (S. pyogenes)	101632
Cas9 (D10A, H840A), NLS, His-Tag (S. pyogenes)	101633
Cas9 (H840A), NLS, His-Tag (S. pyogenes)	101634
Cas9, His-tag (S. pyogenes)	100206
Cas9, NLS, His-Tag (S. pyogenes)	101635
Csm6, His-Tag (E. italicus)	101628
Csm6, His-Tag (T. thermophilus)	101636

Vectors	Catalog#
PD-1 sgRNA-MS2 for CRISPRa (Plasmid)	78091

**Distributed by:**

**CliniSciences Group**

# CliniSciences Group

## Austria

Company: CliniSciences GmbH  
Address: Sternwartestrasse 76, A-1180  
Wien - Austria  
Telephone: +43 720 115 580  
Fax: +43 720 115 577  
Email: [oesterreich@clinisciences.com](mailto:oesterreich@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Belgium

Company: CliniSciences S.R.L  
Address: Avenue Stalingrad 52, 1000  
Brussels - Belgium  
Telephone: +32 2 31 50 800  
Fax: +32 2 31 50 801  
Email: [belgium@clinisciences.com](mailto:belgium@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Denmark

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [danmark@clinisciences.com](mailto:danmark@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Finland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [suomi@clinisciences.com](mailto:suomi@clinisciences.com)  
Web: <https://www.clinisciences.com>



## France

Company: CliniSciences S.A.S  
Address: 74 Rue des Suisses, 92000  
Nanterre- France  
Telephone: +33 9 77 40 09 09  
Fax: +33 9 77 40 10 11  
Email: [info@clinisciences.com](mailto:info@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Germany

Company: Biotrend Chemikalien GmbH  
Address: Wilhelm-Mauser-Str. 41-43,  
50827 Köln - Germany  
Telephone: +49 221 9498 320  
Fax: +49 221 9498 325  
Email: [info@biotrend.com](mailto:info@biotrend.com)  
Web: <https://www.biotrend.com>



## Iceland

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [island@clinisciences.com](mailto:island@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Ireland

Company: CliniSciences Limited  
Address: Ground Floor, 71 lower Baggot street  
Dublin D02 P593 - Ireland  
Telephone: +353 1 6971 146  
Fax: +353 1 6971 147  
Email: [ireland@clinisciences.com](mailto:ireland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Italy

Company: CliniSciences S.r.l  
Address: Via Maremmana inferiore 378  
Roma 00012 Guidonia Montecelio - Italy  
Telephone: +39 06 94 80 56 71  
Fax: +39 06 94 80 00 21  
Email: [italia@clinisciences.com](mailto:italia@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Netherlands

Company: CliniSciences B.V.  
Address: Kraaijenhoffstraat 137A,  
1018RG Amsterdam, Netherlands  
Telephone: +31 85 2082 351  
Fax: +31 85 2082 353  
Email: [nederland@clinisciences.com](mailto:nederland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Norway

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [norge@clinisciences.com](mailto:norge@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Poland

Company: CliniSciences sp.Z.o.o.  
Address: ul. Rotmistrza Witolda Pileckiego 67  
lok. 200 - 02-781 Warszawa -Poland  
Telephone: +48 22 307 0535  
Fax: +48 22 307 0532  
Email: [polska@clinisciences.com](mailto:polska@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Portugal

Company: Quimigen Unipessoal LDA  
Address: Rua Almada Negreiros, Lote 5, Loja 14,  
2615-275 Alverca Do Ribatejo - Portugal  
Telephone: +351 30 8808 050  
Fax: +351 30 8808 052  
Email: [info@quimigen.com](mailto:info@quimigen.com)  
Web: <https://www.quimigen.pt>



## Spain

Company: CliniSciences Lab Solutions  
Address: C/ Hermanos del Moral 13  
(Bajo E), 28019, Madrid - Spain  
Telephone: +34 91 269 40 65  
Fax: +34 91 269 40 74  
Email: [espana@clinisciences.com](mailto:espana@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Sweden

Company: CliniSciences ApS  
Address: Oesterbrogade 226, st. 1,  
Copenhagen, 2100 - Denmark  
Telephone: +45 89 888 349  
Fax: +45 89 884 064  
Email: [sverige@clinisciences.com](mailto:sverige@clinisciences.com)  
Web: <https://www.clinisciences.com>



## Switzerland

Company: CliniSciences Limited  
Address: Marktgasse 18 8302 Kloten -  
Switzerland  
Telephone: +41 (044) 805 76 81  
Fax: +41 (044) 805 76 75  
Email: [switzerland@clinisciences.com](mailto:switzerland@clinisciences.com)  
Web: <https://www.clinisciences.com>



## UK

Company: CliniSciences Limited  
Address: 11 Progress Business center, Whittle  
Parkway, SL1 6DQ Slough- United Kingdom  
Telephone: +44 (0)1753 866 511  
or +44 (0) 330 684 0982  
Fax: +44 (0)1753 208 899  
Email: [uk@clinisciences.com](mailto:uk@clinisciences.com)  
Web: <https://www.clinisciences.com>



## USA

Company: Biotrend Chemicals LLC  
Address: c/o Carr Riggs Ingram,  
500 Grand Boulevard, Suite 210 Miramar  
Beach, FL 32550- USA  
Telephone: +1 850 650 7790  
Fax: +1 850 650 4383  
Email: [info@biotrend-usa.com](mailto:info@biotrend-usa.com)  
Web: <https://www.biotrend-usa.com>

