


Zinc-Stamp™



80849



Zinc-Stamp™ is a chemically labeling reagent that is activated by mobile zinc ion (Zn^{2+}) and thereby covalently attaches a fluorescein tag to proximal proteins. After treatment of cells with Zinc-Stamp™ followed by cell fixation, existence and localization of Zn^{2+} can be analyzed by fluorescent microscopy. In addition, Zinc-Stamp™ enables to analyze and identify proteins involved in Zn^{2+} dynamics and transport.



Features

- Specific for Zn^{2+} . Not reactive with Mn^{2+} , Fe^{2+} , Fe^{3+} , Na^+ , K^+ , Ca^{2+} , and Mg^{2+}
- High S/N ratio
- Low cytotoxicity
- High cell membrane permeability - Just add Zinc-Stamp™ to culture medium
- Labeled proteins can be purified by anti-fluorescein antibody
- Various applications are available; Immunostaining, SDS-PAGE, WB, Proteomics analysis by MS

Product Name	Product Code	Size
Zinc-Stamp™, Mobile Zinc-Responsive Protein Labeling Reagent	FDV-0013A	25μg
	FDV-0013B	3 x 25μg


Distributed by:

CliniSciences Group

Application

Application 1: Imaging of dynamics and localization of mobile Zn²⁺ in cells

The localization of proteins labeled with Zinc-Stamp™ can be visualized by fixing cells. Co-staining with antibodies or organelle-markers can be performed, which enable to analyze the relationships between mobile Zn²⁺ dynamics and protein of interest, such as Zn²⁺ transport proteins.




Cells were treated with S-nitrosocysteine (SNOC; NO-generating reagent) and labeled by Zinc-Stamp™. Cells were fixed by methanol and analyzed with confocal fluorescent microscopy. The responses were observed throughout the cell after 10 minutes of oxidative stress, whereas vesicle-like fluorescent signals were observed after 2 hours of oxidative stress. This indicated that mobile Zn²⁺ was transported to vesicles to suppress toxicity.

Application 3: Analysis of labeled proteins by WB & MS


C6 glioma cells were treated with 200 µM SNOC for 0, 10, and 180 min and incubated with 10 µM Zinc-Stamp™ for 10 min. After lysis of cells, labeled proteins were immunopurified by fluorescein antibody (anti-FL). Analysis were done in Western Blot (A) and Mass Spectrometer (B).

(A) Analysis of individual proteins with western blot



Isolated proteins by anti-FL were separated by SDS-PAGE and detected by western blotting with antibodies of interest (IP-blot). In this case, mitochondrial protein (citrate synthetase), ER protein (calnexin) and Golgi apparatus protein (GRP94) were analyzed. These results indicate mobile Zn²⁺ was released from metallothionein by NO treatment, transiently accumulated in mitochondria within 10 min and finally transported to ER and Golgi apparatus.

Application 2: Analysis of total labeled proteins



C6 glioma cells were treated with 200 µM SNOC, a NO generator, for 0, 15 and 180 min. After NO-stimulation, cells were incubated with 1 µM Zinc-Stamp™ for 10 min in the absence or presence of 400 µM TPEN, a potent Zn-chelator. Cells were lysed by SDS-sample buffer and proteins were separated in SDS-PAGE and detected by fluorescent imager for labeled proteins or CBB staining for total proteins. Band pattern of labeled proteins was dramatically changed by SNOC treatment and fluorescent signals were suppressed in the presence of TPEN.

Your Local Distributor

Funakoshi Co., Ltd. (JAPAN)

Distributed by:

CliniSciences Group

CliniSciences Group

Austria



Company: CliniSciences GmbH
Address: Sternwartestrasse 76, A-1180
Wien - Austria
Telephone: +43 720 115 580
Fax: +43 720 115 577
Email: oesterreich@clinisciences.com
Web: <https://www.clinisciences.com>

Belgium



Company: CliniSciences S.R.L
Address: Avenue Stalingrad 52, 1000
Brussels - Belgium
Telephone: +32 2 31 50 800
Fax: +32 2 31 50 801
Email: belgium@clinisciences.com
Web: <https://www.clinisciences.com>

Denmark



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: denmark@clinisciences.com
Web: <https://www.clinisciences.com>

Finland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: suomi@clinisciences.com
Web: <https://www.clinisciences.com>

France



Company: CliniSciences S.A.S
Address: 74 Rue des Suisses, 92000
Nanterre- France
Telephone: +33 9 77 40 09 09
Fax: +33 9 77 40 10 11
Email: info@clinisciences.com
Web: <https://www.clinisciences.com>

Germany



Company: Biotrend Chemikalien GmbH
Address: Wilhelm-Mauser-Str. 41-43,
50827 Köln - Germany
Telephone: +49 221 9498 320
Fax: +49 221 9498 325
Email: info@biotrend.com
Web: <https://www.biotrend.com>

Iceland



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: island@clinisciences.com
Web: <https://www.clinisciences.com>

Ireland



Company: CliniSciences Limited
Address: Ground Floor, 71 lower Baggot street
Dublin D02 P593 - Ireland
Telephone: +353 1 6971 146
Fax: +353 1 6971 147
Email: ireland@clinisciences.com
Web: <https://www.clinisciences.com>

Italy



Company: CliniSciences S.r.l
Address: Via Maremmana inferiore 378
Roma 00012 Guidonia Montecelio - Italy
Telephone: +39 06 94 80 56 71
Fax: +39 06 94 80 00 21
Email: italia@clinisciences.com
Web: <https://www.clinisciences.com>

Netherlands



Company: CliniSciences B.V.
Address: Kraijenhoffstraat 137A,
1018RG Amsterdam, - Netherlands
Telephone: +31 85 2082 351
Fax: +31 85 2082 353
Email: nederland@clinisciences.com
Web: <https://www.clinisciences.com>

Norway



Company: CliniSciences AS
Address: c/o MerVerdi Munkerudtunet 10
1164 Oslo - Norway
Telephone: +47 21 988 882
Email: norge@clinisciences.com
Web: <https://www.clinisciences.com>

Poland



Company: CliniSciences sp.z.o.o.
Address: ul. Rotmistrza Witolda Pileckiego 67
lok. 200 - 02-781 Warszawa -Poland
Telephone: +48 22 307 0535
Fax: +48 22 307 0532
Email: polska@clinisciences.com
Web: <https://www.clinisciences.com>

Portugal



Company: Quimigen Unipessoal LDA
Address: Rua Almada Negreiros, Lote 5, Loja 14,
2615-275 Alverca Do Ribatejo - Portugal
Telephone: +351 30 8808 050
Fax: +351 30 8808 052
Email: info@quimigen.com
Web: <https://www.quimigen.pt>

Spain



Company: CliniSciences Lab Solutions
Address: C/ Hermanos del Moral 13
(Bajo E), 28019, Madrid - Spain
Telephone: +34 916 750 700
Fax: +34 91 269 40 74
Email: espana@clinisciences.com
Web: <https://www.clinisciences.com>

Sweden



Company: CliniSciences ApS
Address: Oesterbrogade 226, st. 1,
Copenhagen, 2100 - Denmark
Telephone: +45 89 888 349
Fax: +45 89 884 064
Email: [sverige@clinisciences.com](mailto: sverige@clinisciences.com)
Web: <https://www.clinisciences.com>

Switzerland



Company: CliniSciences AG
Address: Fracht Ost Flughafen Kloten
CH-8058 Zürich - Switzerland
Telephone: +41 (044) 805 76 81
Fax: +41 (044) 805 76 75
Email: switzerland@clinisciences.com
Web: <https://www.clinisciences.com>

UK



Company: CliniSciences Limited
Address: 11 Progress Business center, Whittle
Parkway, SL1 6DQ Slough- United Kingdom
Telephone: +44 (0)1753 866 511
or +44 (0) 330 684 0982
Fax: +44 (0)1753 208 899
Email: uk@clinisciences.com
Web: <https://www.clinisciences.com>

USA



Company: Biotrend Chemicals LLC
Address: c/o Carr Riggs Ingram,
500 Grand Boulevard, Suite 210 Miramar
Beach, FL 32550- USA
Telephone: +1 850 650 7790
Fax: +1 850 650 4383
Email: info@biotrend-usa.com
Web: <https://www.biotrend-usa.com>