

**NEOFBS, Fetal Bovine Serum, Collected in South America.**

**Cat# NB-580004A    size: 5 e: 500ml**

**Cat# NB-58-0004B    Size: 100ml**

**General Information**

For decades, fetal bovine serum (FBS) has been an established additive in cell culture media. Its unique growth-promoting properties make it a valuable and essential ingredient for the cultivation of a wide variety of cell types. FBS is a complex mixture of proteins, growth factors, and many other components. It has not yet been possible to produce a synthetic growth additive with comparable performance.

Growing demand of serum with limited availability poses numerous animal welfare issues. With an increasing awareness on this problem, many researchers state that they would like to use less serum – without sacrificing performance.

In collaboration with the University of Applied Science Frankfurt, Neo Biotech developed an alternative to traditional FBS that is more controlled, defined, and less susceptible to the issues of serum supply.

NEOFBS incorporates the key components for optimal cell growth making it possible to significantly reduce the serum content. The combination of these identified components and a lower serum concentration shows a positive influence on cell growth, morphology, and furthermore decreases cell stress while batch-to-batch variations are reduced to a minimum. NEOFBS is supplemented to cell culture media in the same concentration as your standard FBS (e.g., 9 g., 90% medium with 10% FBS).

Product	Origin	Volume	Cat. No.
NEOFBS, Fetal Bovine Serum	Collected in South America	100 ml	<b>NB-58-0004B</b>
		500 ml	<b>NB-58-0004A</b>

**Product Specifications**

Appearance	Clear amber liquid
Storage and shelf life	Store at $\leq -15$ °C. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened, store at 4 °C and use within 4 – 6 weeks.
Shipping conditions	Do not heat inactivate!
Thawing	Frozen (dry ice) Overnight at +2 °C to +8 °C. Swirl gently to homogenize.
	20230731

Ne  
o Biotech  
74 rue des Suiss  
es –  
92000 Nanterr  
e

## Quality Control

The quality of our products starts with the sourcing of raw material, which is only purchased from qualified suppliers and subjected to strict quality controls. After manufacture, the product undergoes extensive testing and only batches that meet all quality criteria are released for sale. Standard tests include pH, osmolality, total protein, albumin, IgG, hemoglobin, endotoxin content, sterility, mycoplasma, virus testing, and cell culture performance.

## Precautions and Disclaimer

This product is for research use only. Please consult the Material Safety Data Sheet (MSDS) for information regarding hazards and safe handling practices.

## Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email ([info@neo-biotech.com](mailto:info@neo-biotech.com)) or phone +33977400909.

For reference only