

# Normal Sheep serum blocking buffer

NB-23-00191-1 (100ml) NB-23-00191-2 (18ml)

### **INTENDED USE:**

Storage: 2-8C

Normal Sheep Serum Blocking Buffer is used in immunohistochemical staining procedures to block non-specific staining when the host species of the secondary antibody is Sheep.

## **COMPONENT:**

Catalog No	Size
NB-23-00191-2	18mL (Ready to use)
NB-23-00191-1	100mL (Ready to use)

Solution contains 0.05% NaN3.

# **RECOMMENDED PROTOCOL:**

- 1. Apply 100uL or enough volume of **Normal Sheep Serum Blocking Buffer** on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
- 2. Blot off Normal Sheep Serum Blocking Buffer Solution (Do not rinse) and add primary antibody

#### **PRECAUTION:**

Wear gloves and take other necessary laboratory safety procedure.

For reference only For Research Use Only.

Neo-Biotech 183, avenue Georges Clémenceau — 92000 Nanterre

