

Normal Bovine serum blocking buffer

NB-23-00180-1 (100ml)
NB-23-00180-2 (18ml)

INTENDED USE:

Storage: 2-8C

Normal Bovine Serum Blocking Buffer is used in immunohistochemical staining procedures to block non-specific staining when the host species of the secondary antibody is bovine (cattle).

COMPONENT:

Catalog No	Size
NB-23-00180-2	18mL (Ready to use)
NB-23-00180-1	100mL (Ready to use)

Solution contains 0.05% NaN₃.

RECOMMENDED PROTOCOL:

1. Apply 100uL or enough volume of **Normal Bovine Serum Blocking Buffer** on tissue before adding primary antibody and incubate for 10 minutes at room temperature.
2. Blot off Normal Bovine Serum Blocking Buffer Solution (Do not rinse) and add primary antibody.

PRECAUTION:

Wear gloves and take other necessary laboratory safety procedure.

*For reference only
For Research Use Only.*

