








Monosan® Tween 20 (50%)

REF / Cat. No.: MON-APP198

Instructions for use

Intended use

50% diluted Tween 20 in water; no added stabilisers; for easier pipetting due to reduced viscosity.

REF	Bestellnummer Catalog Number Reference du catalogue	LOT	Chargenbezeichnung Batch Code Code du lot	 Reizend Irritant Irritant
	Gesundheitsschädlich Harmful Nocif		Giftig Toxic Toxique	 Hersteller / Manufacturer / Fabricant Monosan® Frontstraat 2c 5405 PB Uden The Netherlands Tel: (+31) 413 251115 Fax: (+31) 413 266605 info@monosan.com www.monosan.com
	Verwendbar bis Use By Utiliser jusque			
	Gebrauchsanweisung beachten Consult Instructions for use Consulter les instructions d'utilisation		Lagerungstemperatur Temperature Limitation Limites de température	

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.

07/21 MONOSAN Cat.# MON-APP198

FOR-044 23-10-2019