

Normal Sheep serum blocking buffer #Cat: NB-23-00191-2 Size: 18mL #Cat: NB-23-00191-1 Size: 100mL

Intended use:

Storage: 2-8C

Normal Sheep Serum Blocking Buffer is used in immunohistochemical staining procedures to block non-

specific staining when the host species of the secondary antibody is Sheep.

Component:

Catalog No	Size
NB-23-00191-2	18mL (Ready to use)
NB-23-00191-1	100mL (Ready to use)

Solution contains 0.05% NaN3.

Recommended protocol:

1. Apply 100uL or enough volume of Normal Sheep Serum Blocking Buffer on tissue before adding primary

antibody and incubate for 10 minutes at room temperature.

2. Blot off Normal Sheep Serum Blocking Buffer Solution (Do not rinse) and add primary antibody

3.

Precaution:

Wear gloves and take other necessary laboratory safety procedure.

For reference only For Research Use Only